

## EFFECTIVENESS OF AN ANTACID

---

### PURDUE UNIVERISTY INSTRUMENT VAN PROJECT

#### Determining the Effectiveness of an Antacid Using a pH Meter (Chemistry I Version) (Revised 5/14/96)

#### Introduction

As you may recall, the stomach contains hydrochloric acid. Upset stomach, indigestion, and the much publicized "acid reflux" all occur when the body produces an excess of this acid which travels to more delicate tissues, causing a burning, painful sensation. An antacid is taken to neutralize this excess acid.

Antacids can contain any one or a combination of the following ingredients: aluminum hydroxide, calcium carbonate (commonly known as chalk), various magnesium compounds and sodium bicarbonate. Although these substances are effective, each one has certain disadvantages. For example, sodium bicarbonate loses its effectiveness quickly, and many people must limit their intake of sodium. Magnesium compounds can cause constipation, aluminum hydroxide can act as a laxative, and calcium carbonate has an unpleasant taste.

There are two other important considerations. First, the rate at which the neutralization reaction occurs is very important. It should be relatively fast. The second consideration is the fact that some of these antacids produce carbon dioxide gas as one of the products of the neutralization reaction. In this case, the reaction must occur quickly enough to relieve pain, but slowly enough to allow the gas to dissipate and not build up in one's system.

#### Purpose

In this experiment, you will be using a pH meter to help you determine which brand of antacid is the best at combating indigestion.

#### Materials

3 - 50 mL beakers	1 - 250 mL beaker
1 - 30 mL beaker	stirring rod
4 different brands of antacids	10 mL graduated cylinder
stirring plate and magnetic stir bar	100 mL graduated cylinder
DI water	50 mL buret
mortar and pestle	safety goggles
10 mL 1 M HCl	10 mL 1 M NaOH
100 mL 0.1 M HCl	100 mL 0.1 M NaOH

## EFFECTIVENESS OF AN ANTACID

---

### PURDUE UNIVERISTY INSTRUMENT VAN PROJECT

#### Procedure

##### **Part A- Determination of pH of an antacid dose**

- 1) Grind one dose of an antacid of your choice (use one tablet) using a mortar and pestle.
- 2) Place this dose in a 50 mL beaker. Repeat this process for two more tablets so that you have three 50 mL beakers each containing a single dose.
- 3) Add 10 mL of distilled water to beaker #1. Add 10 mL of 1 M HCl to beaker #2, and add 10 mL of 1 M NaOH to beaker #3.
- 4) Draw or describe what you see occurring in each beaker on your data sheet.
- 5) After the reactions have ceased, check the pH of each solution using a pH meter. Be sure to rinse the electrode with deionized water and blot it dry before transferring the probe to the next beaker.
- 6) Record your data on your data sheet and answer Question #1.

##### **Part B- Determination of the neutralization power of antacids**

- 1) Crush the recommended dose (use 2 tablets) of your antacid using a mortar and pestle. Make the powder as fine as possible, so that it resembles flour.
- 2) Transfer the powder to a 250 mL beaker. Measure 20.0 mL of deionized water and use some of it to rinse the mortar. Pour this into the beaker, along with the remaining water.
- 3) Use a stirring plate to stir the mixture thoroughly until the powder is completely dissolved.
- 4) Describe what happens on your observation sheet.
- 5) Test the pH of the solution using a pH meter and record this value as your initial pH.

## EFFECTIVENESS OF AN ANTACID

---

### PURDUE UNIVERISTY INSTRUMENT VAN PROJECT

- 6) Transfer 5 mL of this solution to a 100 mL beaker.
- 7) Fill a buret with 0.1 M HCl and record the initial volume of the solution. (Remember -- estimate to two decimal places!)
- 8) Place the electrode into the beaker and press the "auto" key to obtain a constant reading of the pH.
- 9) Add one mL of the HCl solution at a time, recording the pH after each addition. Be sure to swirl the beaker carefully after each addition to provide good mixing.
- 10) When the pH reaches 2.0 and remains at or above this level for 30 seconds, stop adding the HCl solution. You have now neutralized all of the antacid and lowered the pH so that it is approximately equal to the pH of the stomach.
- 11) Record the final volume of the HCl solution on your data sheet. This represents the amount of stomach acid a dose of antacid can neutralize.
- 12) Compare your results with those of the other groups and complete your data table.

### Part C- "Cent" and "Scent" analysis

- 1) Record the posted data regarding the costs of various bottles of antacids, and the number of doses contained in each bottle in the Data Table. Calculate the cost per dose for each antacid.
- 2) Record the scents of each brand, noting which is most appealing to you.
- 3) Finally, record the active ingredients for each brand of antacid.

### Questions (to be submitted with your Data Table)

- 1) Which solution from Part A most closely resembles the contents of the stomach?
- 2) What is the relationship between the amount of HCl added and the ability of the antacid to neutralize HCl?

## EFFECTIVENESS OF AN ANTACID

---

### PURDUE UNIVERISTY INSTRUMENT VAN PROJECT

Name: \_\_\_\_\_

Hour: \_\_\_\_\_

#### Antacid Determination Experiment

##### Data: Part A

Antacid Brand Name: \_\_\_\_\_

Beaker #	Observations	pH of solution
1		
2		
3		

Question 1:

##### Data: Part B

Antacid Brand	Observations	Init. Volume HCl (mL)	Final Volume HCl (mL)	Amount of HCl Neutralized (mL)
1.				
2.				
3.				
4.				

## EFFECTIVENESS OF AN ANTACID

---

### PURDUE UNIVERISTY INSTRUMENT VAN PROJECT

#### Data: Part C

Antacid Brand	Bottle Cost	# Doses	Cost per Dose	Scent	Active Ingredients
1.					
2.					
3.					
4.					

Question 2:

#### Conclusion

In addition to submitting your Data Sheet, you are to use your data to analyze which antacid product is the best and most effective. Use the following factors in your analysis: cost per dose, effectiveness, scent and active ingredients. Write this as an article that the general public (parents or other relatives) would be able to understand and use. Be sure to back up your claims with your data! The article does not need to be long, just report your findings and explain.